

# Gradian ICV - Packing List

Item	Part Number	Qty	Packed?
Base Machine   SCXUPT1434			<input type="checkbox"/>
ICV	SCXUPT1413	1	
ICV User Guide (EN)	SCXUPT1097	1	
ICV Quick Start Guide (EN)	SCXUPT0529	1	
ICV Internal Battery	SCXUPT0581	1	
Trolley and Support Arm   SCXUPT1405			<input type="checkbox"/>
ICV Trolley	SCXUPT0463	1	
Circuit Support Arm	SCXUPT0464	1	
Gradian ICV Accessories - Complete Kit   SCXUPT1439			<input type="checkbox"/>
Gradian Accessories 1   SCXUPT1402			
Oxygen Fuel Cell (MOX-3)	SCXUPT1129	1	
Reusable Flow Sensor with Connector (Adult/Peds)	SCXUPT0684	2	
Test Lung (Adult)	SCXUPT0795	1	
Ventilator Bacterial Filter (pack of 2)	SCXUPT0690	2	
Air Inlet Filter	SCXUPT0760	1	
Low Pressure Oxygen Quick-Connect	SCXUPT0691	1	
240 V Humidifier   SCXUPT1164†			
SH530 Servo Controlled Humidifier (AIE) - 240V†	SCXUPT1165	1	
Reusable Humidifier Chamber (Adult/Peds)	SCXUPT0805	1	
Temperature Probe	SCXUPT1179	1	
Power Cord UK - BS 1363 to IEC 60320 C13	SCXUPT0625	2	
Power Cord EU - CEE 7/4 to IEC 60320 C13	SCXUPT1167	2	
Gradian Accessories 2   SCXUPT1403			
Reusable Breathing Circuit (Adult)	SCXUPT0814	2	
Reusable Breathing Circuit (Peds)	SCXUPT0816	2	
NIV Mask (Adult)	SCXUPT0845	1	
HFO Nasal Cannula (Adult)	SCXUPT0849	2	
HMEF (Adult)*	SCXUPT0781	10	
Oxygen Connection Kit   SCXUPT1404			
O2 Hose - DISS(F)-DISS(F)	SCXUPT1102	1	
O2 Hose Adapter - DISS(M)-BS(Probe)	SCXUPT1380	1	
O2 Hose Adapter - DISS(M)-AFNOR(Probe)	SCXUPT1381	1	

\* Consumables Starter Supply

†120V Humidifier Option Available

# Gradian ICV - Optional Accessories

## ICV-K-E | EtCO2 Package | SCXUPT1431 (Optional)

Item	Part Number	Qty	Packed?
ETCO2 Kit			<input type="checkbox"/> (Optional)
Mainstream EtCO2 Module	SCXUPT0582	1	
EtCO2 Airway Adapter (Adult/Peds)*	SCXUPT0682	2	

## ICV-K-N | NICU Package | SCXUPT1438 (Optional)

Item	Part Number	Qty	Packed?
NICU Accessories 1   SCXUPT1429			<input type="checkbox"/> (Optional)
EtCO2 Airway Adapter (Neonatal)*	SCXUPT0697	2	
Reusable Humidifier Chamber (Neonatal)	SCXUPT1103	1	
Reusable Flow Sensor (Neonatal)	SCXUPT0696	2	
NICU Accessories 2   SCXUPT1430			
Reusable Breathing Circuit (Neonatal)	SCXUPT1104	2	
Test Lung (Neonatal)	SCXUPT0698	1	
HFO Nasal Cannula (Neonatal)	SCXUPT0700	2	
HMEF (Neonatal)	SCXUPT0699	10	